

DNA Isolation Protocol

- 1) Get microcentrifuge tube and weigh it
- 2) Allow sample to thaw on ice
- 3) If the sample contains a high moisture content, it will be better to remove the water to obtain a higher DNA yield
- 4) Weigh about 0.5 g of sample in a microcentrifuge tube
- 5) Centrifuge at 10,000 g's for 1 minute
- 6) Remove supernatant
- 7) Reweigh microcentrifuge tube, obtain weight of pellet, ideally 0.25-0.3 grams.
- 8) Repeat steps 3-6 until sample pellet weight is 0.25-0.3 grams
- 9) Break up pellet at bottom of tube using a pipette tip or spatula
- 10) Add ALL contents of Bead Solution Tube to Microcentrifuge tube containing the sample
- 11) Gently vortex to mix.
- 12) **Check Solution S1.** If Solution S1 is precipitated, heat solution to 60°C until dissolved before use.
- 13) Add 60µl of Solution S1 and invert several times or vortex briefly.
- 14) Incubate at 70° C for 10 minutes
- 15) Add 200µl of Solution IRS (Inhibitor Removal Solution). Only required if DNA is to be used for PCR.
- 16) Secure bead tubes horizontally using the Mo Bio Vortex Adapter tube holder for the vortex (cat.13000-V1. Call 1-800-606-6246 for information) or secure tubes horizontally on a flat-bed vortex pad with tape. Vortex at maximum speed for 10 minutes. (See alternative lysis method for less DNA shearing).

- 17) Make sure the 2ml tubes rotate freely in your centrifuge without rubbing. Centrifuge tubes at 10,000 x g for 30 seconds. **CAUTION:** Be sure not to exceed 10,000 x g or tubes may break.
- 18) Transfer the supernatant to a clean microcentrifuge tube (provided).
- 19) **Note:** With 0.25gm of soil and depending upon soil type, expect between 400 to 450µl of supernatant. Supernatant may still contain some soil particles. Discard tube with soil pellet.
- 20) Add 250µl of Solution S2 to the supernatant and vortex for 5 sec. Incubate 4°C for 5 min.
- 21) Centrifuge the tube for 1 minute at 10,000 x g
- 22) Avoiding the pellet, transfer entire volume of supernatant to a clean microcentrifuge tube (provided).
- 23) Add 1.3 ml of Solution S3 to the supernatant (careful, volume touches rim of tube) and vortex for 5 seconds.
- 24) Load approximately 700µl onto a spin filter and centrifuge at 10,000 x g for 1 minute. Discard the flow through, add the remaining supernatant to the spin filter, and centrifuge at 10,000 x g for 1 minute. Repeat until all supernatant has passed through the spin filter. **Note:** A total of three loads for each sample processed is required.
- 25) Add 300µl of Solution S4 (ethanol) and centrifuge for 30 seconds at 10,000 x g
- 26) Discard the flow through.
- 27) Centrifuge again for 1 minute
- 28) Repeat steps 24-26 two times
- 29) Carefully place spin filter in a new clean tube (provided). Avoid splashing any Solution S4 onto the spin filter.
- 30) Add 50µl of Solution S5 to the center of the white filter membrane.
- 31) Incubate at room temperature for 2 minutes

32) Centrifuge for 30 seconds

33) Discard the spin filter. DNA in the tube is now application ready. No further steps are required. We recommend storing DNA frozen (-20°C). Solution S5 contains no EDTA