

DNA Protocol
Gram (-): Negative

- 1.) **Harvest Cells** (max 2×10^9 cells) and **centrifuge** for **10 minutes** at 5000 x g (or 7500 rpm) **remove supernatant**.
- 2.) Added **180 µl** of enzymatic **ATL buffer**.
- 3.) Add **20 µl** of **proteinase K** and **incubate** for **55°C** until cells completely lysis (**30 minutes**) (**Vortex** lightly between incubation to make sure cells are dispersed.)
- 4.) After incubation **vortex** for **15 seconds** and add **200 µl** of **AL buffer** to the sample, mix thoroughly by **vortex** again and incubate for **70°C** for **10 minutes**
- 5.) Add **200 µl** of **Ethanol** (96-100%) to the sample, and (**Vortex**)
- 6.) Pipet mixture for step 6 into a DNeasy spin column placed in a **2ml collection tube**. **Centrifuge** at $>6,000 \times g$ (or 8000 rpm) for **1 minute**. *Discard flow-through and collection tube.*
- 7.) Place DNeasy spin column into a **new 2ml collection tube**. Add **500 µl** of **Buffer AW1** and **centrifuge** for **1 minute** at $> 6000 \times g$ (or 8000 rpm) *Discard flow-through and collection tube.*
- 8.) Place DNeasy spin column into a **new 2ml collection tube**, add **500 µl** of **Buffer AW2**, and **centrifuge** for **3 minutes** at $20,000 \times g$ (or 14,000 rpm) to dry the DNeasy membrane. *Discard flow-through and collection tube.*

- A.) Place DNeasy spin column into a **new 1.5ml microcentrifuge tube**. Pipet **200 µl** of **buffer AE** directly onto the DNeasy membrane. **Incubate** at **room temperature** for **1 minute** and then **centrifuge** for **1 minute** at $>6000 \times g$ (or 8000 rpm) to elute. Save **microcentrifuge tube 1**
- B.) Using a **new microcentrifuge** tube Repeat step 9A. Save **microcentrifuge tube 2**